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



Universal Transfer

**Anvisa Registration:** «80108910051»

**UNIVERSAL TRANSFER** is intended for expert procedures, which must be performed by qualified professionals. The use of the product and surgical techniques are inherent to the professional's training. The use of the product must be performed in a surgical environment and in proper conditions for the health and safety of the patient.

## PRODUCT DESCRIPTION

The Universal Transfer consists of a conical or cylindrical abutment made of plastic material.

## INDICATIONS OF USE

The Universal Transfer is indicated to copy and transfer the position of the component fixed on the implant to a work model where the prosthesis will be made. It aims to transfer the position of the component fixed on the implant in the mouth to a mold (negative) which will then be transformed into a working model (positive).

## OPERATION PRINCIPLE

The working principle of the Universal Transfer is based on the transfer of the position of the component to a mold to prepare the prosthesis.

## HOW TO USE

Stage 01: Definition of the implant feature;  
Stage 02: Definition of the abutment to be used;  
Stage 03: Definition of accessory type to be used;  
Stage 04: Transfer from the implant to a model through the transfer;  
Stage 05: Finalization of the prosthesis on the analogue installed in the model;  
Stage 06: Cemented prosthesis fixation.

## CONTRAINDICATIONS

As long as the material is used properly, there is no contraindication of use.

## SIDE EFFECTS

If the technique is not adequate and the patient is not subjected to the indicated tests, the final result of application of the UNIVERSAL TRANSFER may not succeed, generating product loss or fracture.

## PRECAUTIONS AND RECOMENDATIONS

For the Universal Transfer placement, it is recommended that the professional has a specialization course in the area and prepare a prosthetic execution plan.

Professionals should sterilize all instruments before use, prepare the patient to minimize the risk of contamination and prevent the product comes into contact with any non-sterile object.

The settlement platform of the Universal Transfer that adapt to the Component shall not be changed at all. The professional should be aware of the force exerted when applying the product so as not to damage it.

## WARNING

Single Use, Medical and Hospital Product – Sterile  
Sterilization Process: Gamma Radiation.  
Non-ensured sterilization if the seal is broken.

## TRACEABILITY

All S.I.N. – Sistema de Implante products have sequential batches that allow traceability, which promotes greater safety for the professional qualified to the procedure. Through this batch number, it is possible to know the entire history of the product from the manufacturing process to the distribution time.

## STORAGE

The UNIVERSAL TRANSFER must be stored in a cool and dry place and away from sunlight.

## TRANSPORTATION

The UNIVERSAL TRANSFER must be transported in room temperature, away from direct sunlight, avoiding places where large variations of temperature and humidity occur. The transportation must be carried out properly to avoid falls and it must be carried out in its original package.

## HANDLING CONDITIONS

The UNIVERSAL TRANSFER is a sterile product that should be handled only in a sterile field by properly trained professionals and in appropriate scrubs at the time of the surgical procedure.

## ADDITIONAL INFORMATION

Single use product. Prohibited Reprocessing.  
Exclusive dental use product.

## DISPOSAL OF MATERIAL

The disposal of materials should comply with local hospital regulations and applicable local laws.

## EXPIRATION DATE

Indicated on the label.