INDICATIONS FOR USE

S.I.N. Dental Implant System is intended for placement in the maxillary or mandibular arch to provide support for single-unit or multi-unit restorations. When a one-stage surgical approach is applied, the implant may be immediately loaded when good primary stability is achieved and the functional load is appropriate.

PRODUCT DESCRIPTION

Abutment Cemented is a straight abutment designed for cement-retained prostheses. It is provided in two platform diameters (3.4 and 3.65 mm) for the HE connection and provided for Unitite Slim implants, Unitite implants, Unitite Compact implants, and Strong SW CM implants.

The Abutment Cemented is gold anodized using standard electrolytic passivation processing.

Abutment Cemented is manufactured from Ti-6Al-4V alloy conforming to ASTM F136.

Each abutment is provided with a Ti-6Al-4V alloy abutment screw for attachment to the dental implant.

Provided STERILE. Sterilized by irradiation.

CONTRAINDICATIONS

S.I.N. Dental Implant System is contraindicated in the following conditions:

- The mandibular or maxillary bone quantity and quality is insufficient to provide initial stability to the implant.
- When the site or systemic conditions show inadequate or poor oral hygiene.

- Acute or chronic periodontal infection.
- Chemical dependence.
- Occlusal parafunction.
- Radiation history to the implant site.
- Inappropriate patient for prolonged or complicated oral surgery.
- Inability to build a functional prosthesis.
- Rehabilitation with dental implants is also contraindicated for children, pregnant women and during breastfeeding.
- In cases of immediate loading, inappropriate primary stability of the implant.

WARNINGS

The surgical technique of dental implant installation is highly specialized and the surgical procedure complex, it is recommended that the professionals be technically qualified so that the application of the S.I.N. implants is safe and efficient.

Product is for professional use only.

Product is sterilized by gamma radiation. Sterility is ensured except in cases where the package has been violated or damaged. Do not use if the package is damaged package or after the expiration date. Single use only. Do not re-sterilize.

The reuse or re-sterilization of this product can cause damage to health.

S.I.N. Dental Implant System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) Environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of S.I.N. Dental Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

PRECAUTIONS

Before implant installation, to obtain a predictable long-term outcome, the professional must submit the patient to a detailed and careful medical history, examination, radiographs, laboratory tests, and study models for appropriate planning.

ADVERSE EFFECTS

Loss of the implant and prosthesis is possible due to a number of reasons, including implant contamination, inappropriate surgical technique, poor bone quality, inappropriate oral hygiene, and parafunctional habits (tooth grinding).

SURGICAL COMPLICATIONS

The implant installation surgical procedure may bring risks during and after the surgery, such as: pain, edema, hemorrhage, dehiscence, paresthesia, and infection.

SHIPMENT AND HANDLING

The S.I.N. implants are sent to professionals duly packaged, sealed and sterilized. Therefore, the package must be opened using sterile technique, and must be handled only with sterilized titanium instruments.

ATTENTION

In order to obtain technical support or additional information material about the product, contact: S.I.N. - Sistema de Implante Nacional S.A. Contact details are provided at the end of these instructions.

INSTRUCTIONS FOR USE

After implant installation (1-stage procedure) or after a period of delayed healing (2-stage procedure):

1. Remove the abutment from its packaging and adapt it to the head of the implant with the aid of digital keys (CDH 1220 or CDH1224), handpiece keys (CTH1220 or CTH1224) or ratchet wrenches (CDHC 20 or CDHC 24).

Thread the component onto the implant until it is fully seated onto the implant platform.
Torque to 32 Ncm.

Cemented Abutment Retaining Screws – to attach prosthesis (restoration) to Cemented Abutment

2. Upon completion of the prosthesis (restoration), fit the prosthesis on the Cemented Abutment, and thread the retaining screw onto the abutment until the screw is fully seated.
3. With the aid of 1.2 torque connection wrench and the Prosthetic Torque Wrench, torque the retaining screw to 20 Ncm.

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