CAD-CAM Transfer are intended for expert procedures, which must be performed by qualified professionals. The use of the product and surgical techniques are inherent to the professional's training. The use of the product must be performed in a surgical environment and in proper conditions for the health and safety of the patient.

PROUDCT DESCRIPTION

The CAD-CAM Transfer consists of a conical or square abutment following the seating platforms according to each implant model or prosthetic intermediate. They are manufactured in PEEK and available in NON-STERIL form. They are compatible with Universal Abutments components.

INDICATIONS OF USE

The CAD-CAM Transfer is indicated to copy and transfer the position of the implant to a virtual work model where the prosthesis will be made.

OPERATION PRINCIPLE

Its working principle is based on the transfer of the implant or installed abutment in the oral cavity through an intraoral scan or an extraoral model.

HOW TO USE

- Stage 01:** Definition of the implant or abutment feature;
- Stage 02:** Definition of the CAD-CAM transfer to be used;
- Stage 03:** Screwing of the transfer over the implant or abutment;
- Stage 04:** Scanning of the intraoral transferor or in the model;
- Stage 05:** Prosthesis design in CAD and drilling or printing for finishing (CAM);
- Stage 06:** Fixation of the prosthesis through screw.

ATTENTION

The CAD-CAM Transfer are intended for expert procedures, which must be performed by qualified professionals in implant dentistry. The product must be used in a surgical environment and in proper conditions for the health and safety of the patient.

PRECAUTIONS

For placing the CAD-CAM Transfer it is recommended that the professional has a specialization course in the area and that he elaborates a prosthetic execution plan.

The professional must sterilize the instruments, prepare the patient to minimize the risk of contamination and prevent the product from having contact with any non-sterile object.

The CAD-CAM Transfer that adapts to the implant or prosthetic intermediate should not be altered in any way. The professional must be aware of the force exerted when applying the product so as not to damage it.

RECOMMENDATIONS

For the placement of CAD-CAM Transfer it is recommended that the professional have a specialization course in the area and prepare a prosthetic execution plan.

The S.I.N. does not recommend implant installation in patients with inadequate oral hygiene, uncooperative and unmotivated patients, drug or alcohol abuse, psychoses, chemical dependence, prolonged functional disorders that resist any drug treatment, xerostomia, low immune system, diseases which require the use of steroids regularly, endocrinological diseases, drug allergy, diabetes mellitus, anticoagulants / bleeding diathesis, bruxism, other parafunctional habits, tobacco abuse, installation in children and pregnant women and during breastfeeding.

CONTRAINDICATIONS

As long as the material is used properly, there is no contraindication to use. Transitional product used only to copy and transfer the implant position.

SIDE EFFECTS

The installation recommendations must be followed for the proper functioning of the product, if not, the final result can be compromised generating, loss or fracture of the part. The product can cause transient side effects due to compression of perimplant tissues such as, slight bleeding, edema, pain, discomfort or even infection in case of breaking aseptic barrier.

WARNING

Compatible only with S.I.N.
The product is for single use and cannot be re-sterilized and/or reused.
The reuse or re sterilization of this product may cause contagious infectious disease, deformation and wear of the product.

TRACEABILITY

All S.I.N. – Implant System products have sequential lots that allow traceability, which promotes greater safety for the professional qualified to the procedure. Through this batch number, it is possible to know the entire history of the product from the manufacturing process to the distribution time.

STORAGE

The CAD-CAM Transfer should be stored in a cool dry place at a maximum temperature of 35°C and protected from direct sunlight in their original unopened packaging and should not be damaged.

HANDLING

The CAD-CAM Transfers are non-sterile products that must be sterilized before use and, once sterilized, should only be handled in a sterile environment by properly dressed professionals and wearing appropriate clothing at the time of the surgical procedure.

DISPOSAL OF MATERIAL

The disposal of materials should comply with local hospital regulations and applicable local laws.

TRANSPORTATION

The CAD-CAM Transfer ponents must be transported adequately to avoid falling and stored under a maximum temperature of 35°C, protected from heat and moisture. Carriage must be carried out in its original packaging.

COMPLEMENTARY INFORMATIONS

Single use product. Reprocessing prohibited. Product for exclusive dental use. In the event of an incident caused by the product, the professional must immediately inform the manufacturer.



STERILIZATION

Product provided non-sterile. It must be sterilized in autoclave before use.

1. The product must be enclosed in a steam sterilizable wrap.
2. Steam sterilize in cycles to 121°C at 1 ATM pressure for 30 minutes or to 134°C at 2 ATM pressure for 20 minutes. Drying time 30 minutes.
3. Always accommodate the product in autoclave over a plane surface and away of device walls.
4. Never stack objects or other products.

RECOMMENDATIONS

- a. Sterilize the products in the same day or one day earlier the procedure.
- b. The chemical sterilization is not recommended once some products may cause damages to the product.
- c. Do not use temperature higher than 60°C to drying process.
- d. Do not use dry heat stoves for sterilization of the prosthetic components from S.I.N. - Implant System.

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OBELIS S.A.